



12319 N Mopac Expy, Bldg C, Suite #300, Austin, Tx 78758-2403 (512) 454-0341

## **How to Care for Your Hearing Aids**

- 1. Excessive heat and cold may damage your hearing aids. Never leave them near the stove, in a sunny window, in the car, or other hot places. Do not store the hearing aids in the bathroom. Do not try to dry out the hearing aid in the oven or microwave.**
- 2. Do not wear the hearing aids when using a hair dryer or powertools. Take the hearing aids out while at the barber shop or hair dresser.**
- 3. Do not apply hair spray or other hair products when wearing hearing aids. It can easily damage the microphone.**
- 4. Be careful when wearing hearing aid outdoors or in extremely wet and rainy weather.**
- 5. Never wear the hearing aid while taking a bath or shower, or while swimming.**
- 6. Avoid wearing the hearing aid during strenuous activity, or when you are perspiring excessively.**
- 7. Always carry spare batteries.**
- 8. Batteries should be kept at room temperature. Do not store them in a refrigerator.**
- 9. Keep your hearing aids away from pets and small children.**
- 10. Do not attempt to open up your hearing aid to the internal components. This will void the warranty.**
- 11. Do not wear the hearing aid during a CT-scan, MRI, or similar testing.**